

## **DIAGNOSTIC EQUIPMENT** (INCLUDING, BUT NOT LIMITED TO):

**PULSE OXIMETERS & ACCESSORIES**

**BLOOD PRESSURE MONITORS BLOOD**

**PRESSURE CUFFS & PARTS**

**STETHOSCOPES**

**DOPPLERS**

**THERMOMETERS**

**EAR & EYE SCOPES**

**EKG ELECTRODES**

**AED'S & CABINETS**

**DEFIBRILLATOR PADS & ACCESSORIES**

**DISPOSABLE CO2 SAMPLING LINES**

# **STRATODYNEINC.COM**



## HANDHELD PULSE OXIMETERS & SENSORS

With a huge variety of Handheld Pulse Oximeters available, Stratodyne is certain to have the specific type of system to fit your needs. These units offer dependable results with the flexibility to perform in virtually any setting. Additional features include: adjustable audible & visual alarms, optional data storage and analysis, disposable or rechargeable batteries, units that require no calibration over their lifetime. Stratodyne also carries a full array of all models and sizes of Pulse Oximeter Sensors for all Pulse Oximeters we offer and for any application.



Masimo RAD-57



Nonin® WristOx2



Medsor Handheld



Nellcor SPO2 Monitor



Assurance® Nasal Alar  
SPO2 Sensor



Oximeter Plus 300K1  
Handheld Pulse Oximeter

Masimo Rad-57 Kit w/ SpO2, Adult Rainbow Sensor, RC-01 Cable, Case	5,675.00 EA
Masimo Rad-57 w/ SpO2 only, RC-01 Cable, Adult Soft Boot Sensor, Carry Case	695.00 EA
Nonin 9843 Model Pulse Oximeter / CO2 Detector (Basic System)	1,205.00 EA
Nonin 9847 Model Pulse Oximeter / CO2 Detector (Basic System w/ Full Alarms)	1,305.00 EA
Nonin Carrying Case	49.50 EA
Nonin Carbon Dioxide Sensor	255.00 EA
Nonin Airway Adapter Tubes	7.65 EA 89.05 (12/BX)
Nonin Model 8500 Pulse Oximeter w/ Adult Sensor	585.00 EA
Nonin PalmSAT® Pulse Oximeter	679.00 EA
Nonin PalmSAT® Charger Stand & NiMH Base	296.95 EA
Nonin PalmSAT® Carry Case (Blue)	52.85 EA

**VISIT [STRATODYNEINC.COM](http://STRATODYNEINC.COM) FOR A COMPLETE LISTING OF AVAILABLE SENSORS FOR ALL PULSE OXIMETER BRANDS**

MedSource Hand Held Pulse Oximeter	349.00 EA
SPECTRO 2™ 10 Pulse Oximeter System	419.09 EA
SPECTRO 2™ 30 Pulse Oximeter System	813.65 EA
SPECTRO 2™ Docking Station w/ AC Power	397.85 EA
SPECTRO 2™ Protective Sleeves (Red, Yellow, Green, Blue)	31.89 EA
Smiths/BCI Handheld Pulse Oximeter, Complete	576.75 EA
Smiths/BCI Oximeter Padded Carry Case	57.75 EA
Smiths/BCI Oximeter Protective Rubber Boot	52.55 EA
Nonin /WristOx2 Pulse Oximeter	889.95 EA
Assurance® Nasal Alar SpO2 Sensor Disposable	23.29 EA 558.96 (24/CS)
Nellcore SpO2 Patient Monitor w/ Durasensor Reusable Sensor	649.95 EA
Nellcore SpO2 Patient Monitor w/ Two Sensors, Adult & Pediatric	689.95 EA
Nellcore SpO2 Patient Monitor Carry Case	65.00 EA
ChoiceMed 300K1 Handheld Oximeter	289.00 EA
Fingertip Sensors for the 300K1 Oximeter	94.50 EA

Nonin's® 9840, 9847, 8500, & PalmSAT® Pulse Oximeters w/ Reusable Sensors



SPECTRO 2™ 10 & 30 Pulse Oximeter Systems with Docking Station & Protective Sleeves



Fingertip Probes Available for the 300K1 Handheld Unit

## FINGERTIP PULSE OXIMETERS & ACCESSORIES



Diagnostix 2100 Model



Adimals 2150 Model



Onyx® Vantage 9590



Digit™ Finger Oximeter



Curaplex Fingertip SpO2



Veridian Pediatric Model

Stratodyne currently offers many different brands and styles of Fingertip style Pulse Oximeters. With many adult, pediatric, and combination units in stock, we are certain to carry the model that best suits your needs! Our selection also includes Medical, Sports, and Aviation units. All of our product assure a quick and accurate reading of SpO2 and Pulse Rate measurements.

Visit **StratodyneInc.com** for a complete listing of all our available Fingertip Pulse Oximeters.

Diagnostix 2100 Fingertip Pulse Oximeter (Black,Only)	100.99 EA
Adimals 2150 Fingertip Oximeter (Yellow)	381.05 EA
Onyx® Vantage 9590 Fingertip Pulse Oximeter (Yellow Only)	83.99 EA
Smiths/BCI Digit™ Fingertip Oximeter (Yellow)	381.05 EA
Digit™ Fingertip Oximeter Carrying Case with Belt Loop-Still Nylon Mesh	17.45 EA
Curaplex Fingertip SpO2 Monitor	55.95 EA
Veridian Pediatric Fingertip Pulse Oximeter	59.95 EA

### ALL FINGERTIP PULSE OXIMETERS COME WITH YOUR CHOICE OF A FREE PROTECTIVE CASE OR SKIN

The Elite Fingertip Pulse Oximeter	59.00 EA
The Premier Fingertip Pulse Oximeter	49.00 EA
C21 Fingertip Pulse Oximeter	59.00 EA
C13 Fingertip Pulse Oximeter	49.00 EA
C5 Pediatric Fingertip Pulse Oximeter	118.00 EA
QuickCheck Pro Fingertip Pulse Oximeter	59.00 EA
QuickCheck Fingertip Pulse Oximeter	49.00 EA
Oxi-Go Fingertip Pulse Oximeter	49.00 EA



## FINGERTIP PULSE OXIMETER ACCESSORIES



Keep your Fingertip Pulse Oximeter from getting damaged with this line of accessories from Stratodyne. Choose from our individual protective skins to the complete protective pouch to keep your Pulse Ox free of wear and tear. The protective pouch also keeps your unit under your control by being secured to a belt loop or worn around the neck by utilizing the included lanyard. Many different color options are available.

Fingertip Pulse Oximeter Case w/ Lanyard (Green, Orange, Blue, Black)	2.00 EA
Fingertip Pulse Oximeter Skin (Blue, Black, Red, Translucent White)	2.00 EA

**YOUR CHOICE OF 1 CASE OR 1 SKIN FREE WITH YOUR PURCHASE OF ANY FINGERTIP PULSE OX!!**

## CAPNOGRAPHY

Stratodyne offers several models and types of Capnometers and Capnographs used for short and long term EtCO2 monitoring for adult, pediatric and infant patients. We offer portable, lightweight, and durable units, perfect for operation in any environment. We also provide Capnograph/Pulse Oximeter combination units that are idea to give the user a clear picture of your patient's respiratory status.

EMMA Mainstream EtCO2 Spot Check (3639) w/ Waveform	2,200.00 EA
EMMA Mainstream EtCO2 Spot Check (9632) w/o Waveform	1,600.00 EA
Capnocheck II w/ Battery & Battery Charger	2,744.69 EA
Capnocheck II Capnograph/Digital Oximeter W/SAC Technology	2,895.00 EA
Nonin Lifesense Capnograph	2,929.99 EA
Nonin RespSense Capnograph	2,449.99 EA
Handheld Capnograph w/ Pulse Oximeter	3,985.95 EA
Handheld Capnograph Carrying Case	199.95 EA



Capnometer MA3639



Capnometer MA9639



Capnocheck® II



Lifesense® Model



RespSense™ Model



Handheld Capnograph



## 740 SELECT™ VITAL SIGNS MONITOR

CASMED 740 SELECT™ Vital Sign monitor offers leading edge technology in a rugged, portable, lightweight unit which allows for monitoring of MAXNIBP blood pressure, SpO2 and temperature for patients from neonates to adult even in most challenging situation such as high motion or low perfusion. Many options available. Includes patient cables, adult blood pressure cuff, and startup supplies. The 740 SELECT™ provides an extraordinary selection of standard features that empower the clinician with a high degree of performance, flexibility and customization to address multiple medical environments, enhance patient safety and quickly respond to demanding clinical requirements.

CAS 740 SELECT BP Monitor w/ MAXNIBP	2,429.25 EA
CAS 740 SELECT BP Monitor w/ MAXNIBP, Nellcor SpO2	3,446.25 EA
CAS 740 SELECT BP Monitor w/ MAXNIBP, Masimo or Nellcor SpO2 , Temp	3,653.85 EA
CAS 740 SELECT BP Monitor w/ MAXNIBP, Masimo or Nellcor SpO2 , Temp, (Sidestream or Mainstream)	5,884.50 EA

## MAXTEC V6 VITAL SIGNS MONITOR

The V6 Vital Signs Monitor features an intuitive interface and a compact, portable design. It allows for the monitoring, displaying, reviewing, and storing of a patient's physiological parameters: SpO2, PR, and EtCO2. This V6 also features a large 7" color anti-glare TFT LCD display, a long-lasting rechargeable lithium-ion battery, and a large data storage providing recording capabilities for hundreds of patients.

V6 Vital Signs Monitor w/ ETCO2	2,759.00 EA
V6 Vital Signs Monitor w/ ETCO2, SPO2 & Printer Vet	3,230.00 EA
V6 Vital Signs Monitor w/ ETCO2, SPO2 & NIBP	3,350.00 EA

**CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR A FULL LIST OF ACCESSORIES & PRICING**

## E-SPHYG 2™

This unique backlit LCD display mimics the appearance of a traditional mercury column, but can be easily read in low light conditions. A second numeric LCD displays pressure, pulse, and a host of other vital information. The E-Sphyg 2™ features auto/manual settings. In manual mode, blood pressure is measured using the traditional, time-tested auscultatory method. Auto mode can be used with small adult, adult, and large adult size cuffs to automatically obtain systolic/diastolic pressure and pulse in about one minute. NICAD rechargeable battery provides approximately 200 measurements on a single charge. AC adapter included. Three year warranty.

Diagnostix E-Sphyg 2™ Multi Cuff Kit, Desk or Wall Model	423.25 EA
Diagnostix E-Sphyg 2™ Multi Cuff Kit, Mobile Model	506.09 EA

## ADVIEW™ BLOOD PRESSURE SYSTEMS & CUFFS

Designed with the idea that no two patients are alike, this is the first clinical grade automated blood pressure device to offer a manual back-up mode. By combining the best of blood pressure research and technology with the capability to apply traditional methods of blood pressure measurement, the ADView uses the most reliable means of capturing blood pressure data.

ADView™ BP Module	597.85 EA
ADView™ BP & Thermometry Modules	1,155.35 EA
ADView™ BP & SPO2 Modules	1,187.45 EA
ADView™ BP, SPO2, & Thermometry Modules	1,804.25 EA

## SPOT VITAL SIGNS™

Welch Allyn's Spot Vital Signs™ device delivers vital signs in seconds. The basic unit offers noninvasive blood pressure, pulse rate and mean arterial pressure (MAP). Spot Vital Signs additional parameters include SureTemp™ thermometry and/or Nellcor® Pulse Oximetry.

Sport Vital Signs™ BIBP Unit Only	1,105.95 EA
Sport Vital Signs™ NIBP & Nellcor Pulse Oximeter Unit	2,165.95 EA
Sport Vital Signs™ NIBP, Nellcor Pulse Oximeter, & Thermometry Unit	2,562.95 EA

**CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR A FULL LIST OF OPTIONS & PRICING**

## ADVANTAGE™ DIGITAL BLOOD PRESSURE MONITORS

With Smart Logic™ Technology and an advanced memory, the entire ADVANTAGE™ product line has become the most popular arm and wrist digital home blood pressure monitors. Varying in styles, each monitor measures systolic, diastolic pressure, and pulse.

Advantage™ 6015N Digital Wrist BP Monitor	46.95 EA
Advantage™ Ultra Wrist Digital BP Monitor 6016N	76.65 EA
Advantage™ 6012N Semi-Automatic BP Monitor	38.25 EA
Advantage™ 6021N Automatic Digital BP Monitor	66.95 EA
Advantage™ Plus 6022N Automatic Digital BP Monitor	76.95 EA



## PROSPHYG™DESIGNER SERIES BLOOD PRESSURE CUFFS

Our ADC Blood Pressure Cuffs feature a black enamel 300mmHG no-pin stop manometer that provides day-in, day-out dependability. ADC's proprietary Size Guide marking system prevents mis-cuffing. Each series includes a chrome-plated brass air release valve. 768 Series Includes a color-matched, dual pocket carry case. Pro's Combo II™ SR Series includes premium nylon carry case with an Adscope™ Sprague 641.




Prosphyg™ 760 Series (Adult, Sm. Adult, Child)	36.29 EA
Prosphyg™ 768 Series (Adult, Sm. Adult, Child)	35.65 EA
Pro's Combo II™ SR (Adult, Sm. Adult, Child)	48.19 EA

**THEMES AVAILABLE ARE: BLUE SWIRLY, MEDICAL, ANIMALS, BREAST CANCER AWARENESS, & PUZZLE PIECES**

## ADC PROSPHYG™760, 775, 770 SERIES

Our Prosphyg™ feature a black enamel 300mmHG no-pin stop manometer that provides day-in, day-out dependability. ADC's proprietary Size Guide marking system prevents mis-cuffing. Each series includes a chrome-plated brass air release or deflation valve and latex inflation bladder. 775& 770 Models include a bladder & bulb, as well as artery label & gauge holder. 770 Model also includes gray cotton cuff.



Prosphyg™ 760 Series (Sm. Adult, Child, Infant) (Black)	29.40 EA
Prosphyg™ 760 Series (Adult) (Black, Royal, Gray, Red, Orange)	29.40 EA
Prosphyg™ 760 Series (Lg. Adult) (Black)	34.99 EA
Prosphyg™ 760 Series (Thigh) (Black)	39.45 EA
Prosphyg™ 775 Series (Child, Sm. Adult, Adult, Lg. Adult) (Navy)	18.95 EA
Prosphyg™ 775 Series (Thigh) (Navy)	24.95 EA
Prosphyg™ 770 Series (Adult) (Gray Cotton)	17.79 EA

## WELCH ALLYN ANEROIDS & CUFF KITS

The Ultimate hand aneroid sphygmomanometers that provide full flexibility of reading positions with ease of use. Precise control for inflation and deflation. The DS66 Model is updated with an ergonomic design and 1 finger operation. Using a trigger style release valve to provide smooth, even, and far superior deflation. Trigger moves to locked position after complete deflation. Multi-Cuff Kit includes an adult, large adult, and child cuff with a dual compartment zippered case. Certified 4 Cuff Kits come with a hand aneroid, a handy two compartment case, and either a adult or pediatric cuff assortment.



DS66 Trigger Aneroid      Certified 4 Cuff Kit      Certified Hand Aneroid

Adult DS66 Trigger Aneroid (Complete)	189.00 EA
DS66 Trigger Multi Cuff Kit	286.79 EA
Certified Adult Hand Aneroid (Complete)	265.35 EA
Adult Certified 4 Cuff Kit	331.95 EA
Pediatric Certified 4 Cuff Kit	296.49 EA

**CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR A COMPLETE OFFERING & PRICING OF REPLACEMENT PARTS AVAILABLE**

## ADC PROSPHYG™DIAGNOSTIX & E-SPHYG™ SERIES

The Diagnostix series are precision crafted 300mmHg no-pin stop German movement tested to 200,000 cycles. Options varying per each Series type include: chrome plated manometers, swivel, or ADCUFF™ cuffing for proper positioning, latex or ADFLOW™ bladder & bulb for easy inflation or deflation. The E-Sphyg™ replaces your mechanical gauge with a virtually indestructible, digital instrument that can stand up to the toughest use, with auto off function, swiveling cuff clip, and pulse rate display.



**Displayed Clockwise From Left:** 703 Series, 700 Series, 750 Series & E-Sphyg™

Diagnostix™ 703 Series (Adult, Sm. Adult, Child)	82.79 EA
Diagnostix™ 700 Series (Adult, Sm. Adult, Child, Infant) (Black)	40.95 EA
Diagnostix™ 700 Series (Lg. Adult) (Black)	46.65 EA
Diagnostix™ 700 Series (Thigh) (Black)	50.25 EA
Diagnostix™ 750 Series Wall Mounted	98.45 EA
Cuff Storage Basket for Diagnostix™ 750 Series	14.39 EA
E-Sphyg™ Complete Sphygmomanometer (Adult) (Black or Navy)	67.25 EA

## BLOOD PRESSURE MEDIC KITS

Featuring a variety of medic blood kits to choose from. Each kit includes a variety of either 3, 4, or 5 standard sized cuffs. Each kit also includes a durable, multi-pocketed case to keep the kit components organized and protected. Curaplex Kit includes a chrome plated aneroid gauge. Multikuf™ Kit includes an ambidextrous palm gauge with trigger valve. System 5™ Kit includes DIAGNOSTIX™ A804 palm style manometer and an ambidextrous trigger release air valve.

Curaplex B.P. Kit (3) (Blue or Orange)	57.50 EA
Curaplex B.P. Kit (5) (Blue or Orange)	67.50 EA
Multikuf™ BP System (3) (Black, Navy, or Orange)	127.45 EA
Multikuf™ BP System (4) (Black, Navy, or Orange)	142.25 EA
System 5™ Multicuff BP Kit (Navy or Orange)	154.69 EA



Curaplex B.P. Kits



ADC Multikuf™



System 5™

## STETHOSCOPES & COMBINATION KITS

Stratodyne carries a large inventory of Stethoscopes as well as Stethoscope Combination Kits. We have many models and styles to choose from including models specifically for pediatric patients. Our stethoscope selection varies from lightweight aluminum to surgical stainless steel. Choose between various tubing lengths as well as bell construction for patient comfort and performance. We also carry a line of Electronic Stethoscopes for improved sound quality. Combination kits include Stethoscope, BP Cuff & Manometer, Penlight, Shears, and carrying case.



Pro's Combo IV Kit



Sterling Series ProKit



Sprague Rappaport Type



Proscopes™ Stethoscope



Pediatric/Infant Adscopes™



Adscope™ Lite 609



3M™ Littmann™ Classic II S.E.



Classic II Pediatric Stethoscope



3M™ Littmann™ Electronic Stethoscope

<b>Pro's Combo II Pocket Aneroid Kit</b> (Black, Orange, Burgundy, Red, Magenta, Teal, Purple, Royal, Navy, Dark Green)	42.75 EA
<b>Pro's Combo IV Fanny Pack Kit</b> (Black, Royal, Navy, Orange)	60.49 EA
<b>Sterling Series ProKit</b> (Black, Navy, Royal, Burgundy, Hunter Green, Magenta, Purple, Teal)	9.95 EA
<b>Sterling Series Sprague Rappaport Stethoscope</b> (Black, Navy, Royal, Hunter Green, Magenta, Orange, Red, Yellow)	9.95 EA
<b>Curaplex Sprague-Rappaport Stethoscope</b> (Black, Blue, Green, Red, Yellow, Orange)	9.95 EA
<b>ADC Sprague-Rappaport Type Stethoscope</b> (Navy, Black, Burgundy, Red, Neon Orange, Neon Yellow, Teal, Royal, Frosted Glacier, Frosted Seafoam, Frosted Lilac Frosted Magenta)	11.68 EA
<b>ADC Adscope® 606 Ultra-lite Cardiology Stethoscope</b> (Black or Tactical)	75.75 EA
<b>ADC Platinum Edition Adscope-Lite™ Stethoscope</b> (Black)	52.35 EA
<b>ADC Adscope™ 600 Platinum Edition Stethoscope</b> (Black)	130.59 EA
<b>ADC Pediatric or Infant Adscopes™</b> (Black, Light Blue)	53.95 EA
<b>ADC Adscope™ Lite 609</b> (Black, Navy, Burgundy)	21.95 EA
<b>ADC Adscope™ 603 Stethoscope</b> (Black, Burgundy, Royal, Red)	36.59 EA
<b>ADC Adscope™ 601 Stethoscope</b> (Black, Navy)	99.75 EA
<b>ADC Proscopes™ Stethoscope</b> (Nurse Scope) (Navy, Black, Burgundy)	5.95 EA
<b>ADC Proscopes™ Stethoscope</b> (Dual Head) (Navy, Black, Red, Royal)	7.79 EA
<b>ADC Proscope™ 675 Dual Head Newborn Scope</b> (Pink)	14.79 EA
<b>ADC Proscope™ 675 Dual Head Pediatric Scope</b> (Pink or Gray)	9.95 EA
<b>ADC Proscope™ Teaching Dual Head Scope</b> (Black)	18.65 EA
<b>3M Littman Cardiology IV</b> (22" or 27") (Black, Burgundy, Navy, Green)	243.65 EA
<b>3M Littman Master Cardiology</b> (Black, Burgundy, Navy, Hunter Green, Plum)	232.75 EA
<b>3M™ Littmann™ Electronic Stethoscope Model 3100</b> (Black, Burgundy, Navy)	454.49 EA
<b>3M Littman Master Classic II</b> (Black, Burgundy, Navy, Caribbean Blue)	130.59 EA
<b>3M Littman Classic II S.E.</b> (Black, Navy, Cell Blue)	82.25 EA
<b>3M Littman Lightweight II S.E.</b> (Black, Burgundy, Caribbean Blue, Lilac, Cell Blue, Seafoam Green)	55.59 EA
<b>3M Littman Select</b> (Black, Burgundy, Caribbean Blue, Royal)	60.79 EA
<b>3M Littman Infant Classic II</b> (Black)	101.39 EA
<b>3M Littman Pediatric Classic II</b> (Black or Navy)	101.39 EA
<b>3M Littman Classic III Stethoscope</b> (Black, Navy, Raspberry, Burgundy, Turquoise, Lemon-Lime)	102.95 EA
<b>S.T.H. Stethoscope Tape Holder</b> (Neon Green, Pink, Red, Royal, Black)	5.95 EA

**CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR ALL AVAILABLE SPECIFICATIONS & PRICING**

## DIGITAL DOPPLERS & DOPPLER SCOPES

The DigiDop & DigiDop II utilize digital signal processing to reduce unwanted background noise and let the real Doppler sounds through bright and clear. The DigiDop II features extended length probes for up to 56% increased sensitivity at shallow depths. The smaller size and shape were designed to comfortably fit in your hand or coat pocket. The CardioBeat Doppler Scope features a patented noise reduction feature to allow even the weakest pulse to resonate clearly. It also features an enhanced detection mode for increased sensitivity while screening more difficult patients.

DigiDop (Non-Display) (Non-Rechargeable) w/ 3MHz or 8MHz Probe	524.95 EA
DigiDop (Non-Display) (Rechargeable) w/ 3MHz or 8MHz Probe	578.95 EA
DigiDop (Backlit Display) (Non-Rechargeable) w/ 3MHz or 8MHz Probe	724.95 EA
DigiDop (Backlit Display) (Rechargeable) w/ 3MHz or 8MHz Probe	772.95 EA
DigiDop II (Display) w/ 3MHz or 8MHz Probe	725.00 EA
DigiDop II (Display) (Recharger) w/ 3MHz or 8MHz Probe	775.00 EA
CardioBeat Doppler Scope	696.95 EA
CardioBeat Doppler Scope Gel (60g)	2.20 EA



DigiDop



DigiDop II



CardioBeat



## THERMOMETERS

Stratodyne provides a full and complete offering of Thermometers for critical and clinical patient measuring applications that require higher precision accuracy. Featuring many different styles and models, we are certain to have the right thermometer for your needs. Oral, Rectal, axillary and multi-probe, as well as ear and forehead scanning is available. Options are available for replacement probe covers or models that are cleaned with alcohol. We also feature 3M's Tempa-Dots, which are a single use clinical thermometer used to provide accurate and reliable measurements. They are sterile and individually wrapped.

**VISIT STRATODYNEINC.COM TODAY FOR A FULL PRODUCT OFFERING!!**

MedSource™ Non-Contact Thermometer	69.95 EA
Pro 6000 Thermoscan Tympanic Thermometer w/ Small Cradle	250.00 EA
Pro 6000 Thermoscan Tympanic Thermometer w/ Large Cradle	255.00 EA
Charging Station for Pro 6000 Thermoscan Tympanic Thermometer	149.00 EA
GENIUS™ 2 Tympanic Thermometer w/ Base	389.95 EA
SureTemp® Plus 690 Thermometer	296.95 EA
SureTemp® Plus 690 Thermometer w/ Holder	310.95 EA
Oral Probe/Wellness Kit or Rectal Probe/Wellness Kit for SureTemp® Plus 690 Thermometer	89.00 EA
SureTemp® Plus 692 Thermometer (Oral/Axillary Probe & 4' Cord)	486.00 EA
Adtemp™ 427 Temple Touch Thermometer	12.95 EA
Veridian 9-Second Digital Thermometer (Rigid or Flexible Tip, Storage Case, 5 Probe Covers)	9.95 EA
Veridian 60-Second Digital Thermometer	3.35 EA
Veridian Touch-Free Infrared Thermometer	49.95 EA
Veridian Non-Contact Infrared Thermometer	74.95 EA
MOBI DualScan Ultra 2 in 1 Digital Thermometer	28.99 EA
ADC Adtemp™ 421 Digital Ear Thermometer	34.79 EA
ADC Hypothermia Thermometer	10.75 EA
ADC Adtemp™ II 413 Digital Thermometer	5.59 EA
3M Tempa-Dots (100 Per Box)	15.99 BX

**PLEASE CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR PRICING  
ON REPLACEMENT PROBE COVERS AND ACCESSORIES**



MedSource™  
Non-Contact



Pro 6000 Thermoscan  
Tympanic



SureTemp® Plus 690



SureTemp® Plus 692



427 Temple Touch



9 Second Digital



Touch Free Infrared



60 Second Digital



Tempa-Dots



Adtemp™ II 413  
Digital



ADC Hypothermia



Adtemp™ 421  
Digital Ear



DualScan Ultra 2 In 1

## EAR & EYE SCOPES

ADC'S Standard, Economy, and Pocket Otoscope/Ophthalmoscope Sets allow for a precise visual inspection of the eyes and inner ears of a patient when required. The **Otoscope** features: obstruction free light transmission; a wide angle, removable 2.5x viewing lens; a sealed system with included insufflator fitting for pneumatic otoscopy, an Otoscope head made from impact resistance polymer, interchangeable battery handle with the Ophthalmoscope, and tubes of (5) 2.5mm and (5) 4.0mm specula. The **Ophthalmoscope** features: (5) aperture selections, (17) lenses from -20 to +20 diopters for better resolution, an illuminated dial, built-in lens dust cover, silicon coated protection of glass components, rubber brow rest to prevent scratching, and a battery handle that interchanges with the Otoscope.

Otoscope/Ophthalmoscope Standard Set	237.35 EA
Economy Pocket Otoscope/Ophthalmoscope Set	205.95 EA
Pocket Otoscope	116.65 EA



Standard Set



Economy Set



Pocket Otoscope

## EKG ELECTRODES

## ECG SNAP ELECTRODES

Dynarex Electrodes feature Ag/ AgCl coating on electrode surface for effective signal transmission.

3p-Pack	0.50 PK	96.25 (200PK/CS)
4-Pack	0.60 PK	87.30 (150PK/CS)
50-Pack	5.35 PK	62.35 (12PK/CS)

## FOAM ELECTRODES

Features a "wet" gel for maximum adhesion & electrical contact.

3p-Pack	0.65 PK	126.10-200PK/CS
4-Pack	0.89 PK	126.10-150PK/CS
30-Pack	10.50 PK	203.70-20PK/CS
10-Pack	2.10 PK	122.25-60PK/CS

## CLEARTRACE™

The original clear tape (Neoderm®) electrode with conductive adhesive gel.

Box (30)	8.98 BX	174.20 (600/CS)
Pouch (30)	12.75 PKG	247.35 (600/CS)

## POSITRACE

Latex free & low profile foam for increase patient comfort.

Box (30)	8.49 BX	164.70 (600/CS)
Pouch (30)	8.19 PKG	158.85 (600/CS)

## PEDIATRIC HUGGABLES

Latex free with perforated tape along with small size for comfort. A gel pocket stabilizes the electrode.

Box (30)	8.75 BX	169.75 (600/CS)
----------	---------	-----------------

## INSTATRACE

Latex free with closed cell foam and a pull tab for easy application. We gel formula provides quick signal retention.

Box (30)	10.85 BX	210.49 (600/CS)
Pouch (30)	9.75 PKG	189.15 (600/CS)

## ULTRATRACE

Features we gel application for strong signals and stable baselines, Neoderm® tape, and easy application tab.

Box (30)	8.05 BX	156.15 (600/CS)
----------	---------	-----------------

## EXACTRACE

Latex free with closed cell foam to keep fluids from penetrating. Wet gel for quick signal retention and a pull tab for easy application.

Box (30)	11.15 BX	216.30 (600/CS)
Pouch (30)	10.69 PKG	207.39 (600/CS)

## MEDI-TRACE™

High performance in stress, Holter, & diaphoretic applications. Conductive adhesive hydrogel enables consistency.

30-Pack	7.05 PK	136.75 (600/CS)
3-Pack	0.90 PK	174.60 (200PK/CS)
5-Pack	1.50 PK	174.60 (120PK/CS)

## BIOTAC ULTRA

Feature long term adhesion, low profile, & resistant closed cell foam. (3) Electrodes per package.

Box (150)	29.95 BX	116.19 (600/CS)
-----------	----------	-----------------

## MEDI-TRACE MINI™

Features a distinct & size and a soft closed cell foam for maximum adhesion. (3) Electrodes per package.

Box (60)	13.28 BX	128.80 (600/CS)
----------	----------	-----------------

## ELECTROBLUE® STRESS

Impervious to fluids, tear resistant, and conformable to body contours. Medical grade adhesive designed for demanding conditions.

3-Pack	30.95 (75/BX)	120.05 (300/CS)
22.50 (50/PK)	218.25 (500/BX)	854.95 (2000/CS)

## EXCEL STRESS FOAM

Radiotranslucent solid gel electrodes that are impervious to fluids, tears, and conform to body contours. Conductive hydrogel provides "artifact free" tracings.

12.49(50/PK)	121.15 (500/BX)	489.60 (2000/CS)
--------------	-----------------	------------------

## MEDITRACE® 230

Features a conductive adhesive hydrogel, a smaller size for convenient lead placement, and high quality foam substrate to resist fluids and ensure trace quality.

30-Pack	5.79 PK	112.35 (20PK/CS)
---------	---------	------------------

## RED DOT SOFT CLOTH

Same size, round-style, soft-cloth, electrodes designed for general or long-term clinical situations.

30-Bag	20.29 BG	393.65 (20 BG/CS)
--------	----------	-------------------

## RED DOT FOAM

3M foam composition. Each electrode includes an abradar pad to ease patient preparation.

3-Pack	31.75 (20PK/BX)	307.95 (10BX/CS)
50-Pack	21.95 PK	425.85 (20PK/CS)

## RED DOT DIAPHORETIC

Diaphoretic Foam and cloth electrodes packaged 50 per each bag. These high quality electrodes provide exceptional performance.

Foam	21.99 (50/PK)	426.60 (20PK/CS)
Cloth	23.45 (50/PK)	454.90 (20PK/CS)

## CARDIOTRACE UNIVERSAL

Low profile foam, solid gel, "snap on" electrodes feature a lift tab for effortless removal from the liner card and from the patient's skin.

10-Pack	2.35 PK	113.95 (50PK/CS)
---------	---------	------------------

## AMBU BLUE SENSOR®

Offset foam composition maintains a steady uninterrupted signal while maintaining patient comfort.

1.95 (LG. 3-PK)	2.35 (LG. 4-PK)	2.65 (LG. 5-PK)
5.85 (LG. 10-PK)	12.10 (LG. 25-PK)	4.99 (MD. 10-PK)
15.75 (MD. 50-PK)	14.59 (SM. 25-BG)	

## CENTER STUD (1.5) ECG

Highly conductive solid gel with strong adhesion to ensure good signal quality during short monitoring applications.

5-Pack	10.85 (12/PK/BX)	105.20 (10/BX/CS)
10-Pack	10.39 (6/PK/BX)	100.59 (10/BX/CS)
50-Pack	9.60 PK	111.79 (12/PK/CS)

## VERMED EMS

Specifically for EMS, insures quick placement & fast, accurate readings. Hypoallergenic & latex-free.

3/Strip 75 Per Box	25.89 BX	200.79 (8BX/CS)
4/Strip 100 Per Box	38.95 BX	226.49 (6BX/CS)
10 Pack For 12 Lead ECG Kit		188.35 (50PK/CS)

## VERMED 4 LEAD EMS

1 1/2" Round with aggressive adhesive to stick on diaphoretic patients.

4/Strip	133.69 (600/CS)
---------	-----------------

## VERMED PEDIATRIC EMS

3 and 4 lead pediatric electrodes with mild adhesive and lower chloride gel for young patients.

4 Per Strip	32.85 (100/BX)	191.25 (600/CS)
-------------	----------------	-----------------

VISIT [STRATODYNEINC.COM](http://STRATODYNEINC.COM) FOR A  
COMPLETE LISTING OF EKG  
ELECTRODES ALONG WITH ALL SIZING  
& QUANTITY OPTIONS!!



## AUTOMATIC EXTERNAL DEFIBRILLATORS

Stratodyne carries a huge variety of AEDs in order to ensure that we have something to meet your specific needs and requirements. With a supply consisting of multiple manufacturers and many models per manufacturer to choose from, we have put together a product offering to make sure we can meet or exceed your expectations. We also carry a wide array of Recertified AEDs and accessories. The table below shows only a portion of what we have to offer, please visit **StratodyneInc.com** or contact a Stratodyne representative for a full product offering, including accessories for each specific system.



Philips HeartStart MRx



Philips HeartStart FR3



Philips HeartStart Onsite



Cardiac Science  
Powerheart G3



Cardiac Science  
Powerheart G3 PRO



Defibtech  
Lifeline



Defibtech  
Lifeline Pro

Philips HeartStart MRx Defibrillator	Contact for Pricing
Philips HeartStart FR3 Defibrillator (Text Bundle)	2,530.00 EA
Philips HeartStart FR3 Defibrillator (ECG Bundle)	2,990.00 EA
Philips HeartStart FRx Defibrillator	1,691.00 EA
Philips HeartStart OnSite Defibrillator	1,290.00 EA
Philips HeartStart OnSite Defibrillator (Redi-Pack)	1,391.99 EA
Cardiac Science Powerheart GC AED Automatic	1,595.00 EA
Cardiac Science Powerheart GC AED Semi-Automatic	1,595.00 EA
Cardiac Science Powerheart G3 Pro AED	2,549.95 EA
Cardiac Science Powerheart G5 AED (Fully or Semi Automatic)	1,695.00 EA
-Available in English Only or English & Spanish	
Defibtech Lifeline AED	1,245.00 EA
Defibtech Lifeline Automatic AED	1,295.00 EA
Defibtech Lifeline View AED	1,595.00 EA
Physio-Control LIFEPAK EXPRESS® AED	1,695.00 EA
Zoll AED Plus Defibrillator	1,699.00 EA


**CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR A FULL LISTING  
OF AVAILABLE ACCESSORIES AND PRICING**



Philips HeartStart FRx



Cardiac Science Powerheart G5  
(Auto & Semi Auto)



Physio-Control  
LIFEPAK EXPRESS



Zoll AED Plus

## DEFIBRILLATOR CABINETS

Our well constructed AED Wall Cabinets provide the best way to store your AED in an easily accessible and highly visible location. These units can be customized with several different options including with or without an audible alarm.

Cardiac Science Surface Mount AED Wall Cabinet w/ Alarm	189.00 EA
Cardiac Science Surface Mount AED Wall Cabinet w/out Alarm	249.00 EA
Zoll AED Standard Metal Wall Cabinet with Alarm	205.00 EA
Philips AED Cabinet with Alarm	228.65 EA
Physio-Control LifePak® Wall Cabinet w/ Alarm (Semi-Recessed)	315.15 EA



Cardiac Science



Zoll



Philips



LifePak®

## AUTOMATIC ELECTRONIC DEFIBRILATORS OPTIONS &amp; ACCESSORIES



Disposable CO2 Lines



OMNISTREAM™ Airway Adaptor




Microstream FilterLine® H Set



Smart SureLine™ Plus



SureStream™ WaterTrap



Smart OmniLine™ Plus

## DISPOSABLE CO2 SAMPLING LINES

Stratodyne offers multiple products to be used for CO2 sampling. The various features of each product include models for Nasal and Oral sampling; adult, intermediate, and pediatric style units; and units specially designed for short or long term use.

Disposable CO2 Sampling Line (Adult, w/ O2 Tubing)	13.29 EA	
Disposable CO2 Sampling Line (Pediatric, No O2 Tubing)	13.79 EA	334.00 (25/BX)
Disposable CO2 Sampling Line (Pediatric, w/ O2 Tubing)	14.69 EA	367.25 (25/BX)
Disposable CO2 Sampling Line (Adult, w/ O2 Connector)	12.59 EA	1,221.25 (100/BX)
Disposable CO2 Sampling Line (Adult, Non-Humidified)	10.25 EA	994.25 (100/BX)
OMNISTREAM™ Airway Adaptor	5.35 EA	133.75 (25/BX)
Microstream FilterLine® H Set (Adult/Pediatric)	16.75 EA	418.75 (25/BX)
Microstream FilterLine® H Set (Infant/Neonate)	17.05 EA	413.45 (25/BX)
Smart CapnoLine® Plus Long (Adult/Intermediate) (w/ O2 Delivery)	25.75 EA	
Microstream FilterLine® SET (Adult/Pediatric)	11.25 EA	
Smart SureLine™ Plus Disposable (Pediatric) (w/ O2 Tubing)	18.49 EA	462.25 (25/BX)
SureStream™ WaterTrap	17.09 EA	427.25 (25/BX)
Smart OmniLine™ (Pediatric) (w/ O2 Tubing)	9.50 EA	237.50 (25/BX)
Smart OmniLine™ (Adult/Intermediate) (w/ O2 Connector)	9.40 EA	235.00 (25/BX)
Smart OmniLine™ (Adult/Intermediate) (w/ O2 Tubing)	9.90 EA	247.50 (25/BX)

## DEFIBRILATOR PADS


Stratodyne provides Defibrillator pads from every major manufacturer. We are certain to carry the specific pads that your operation requires. We also carry all major brands of Defibrillator accessories. This includes: Recording Paper, Replacement Batteries, Paddle Gel, Skin Preparation Equipment, and Electrode Cables. Please visit [StratodyneInc.com](http://StratodyneInc.com) for a full product listing complete with specifications and additional options



ZOLL Stat Padz II



ZOLL CPR-D Padz

Physio-Control LifePak®  
Edge Quick Combo  
Multifunction PadsPhysio-Control Quik-  
Combo Defib Pads

PadPro™



Kendall

Cardiac Science  
Adult & PediatricHeartstart MRx  
Multifunction

ZOLL CPR-D Padz Kit	170.95 EA	
ZOLL Stat Padz II (Multifunction)	59.95 PR	
ZOLL CPR Stat Padz (HVP Multi-Function) (M, R, & E Series)	75.95 PR	
ZOLL Pro Padz Radiolucent Solid Gel Electrodes	66.95 PR	
ZOLL Pedi Padz Multi-Function Solid Gel Electrodes	60.95 PR	
ZOLL Pedi Padz II Multi-Function Electrodes for AED Plus & Pro	96.95 PR	
Physio-Control LifePak® Edge Quick Combo Multifunction (RTS)	Adult: 49.95 PR	Ped.: 45.95 PR
Physio-Control LifePak® Edge Fast Patch Plus Defib Pads	27.95 PR	
Physio-Control Quik-Combo Defib Pads (Redi-pack)	46.99 PR	
Physio-Control Infant/Child Reduced Energy Pads	115.50 PR	
PadPro™ Defib Pads (Zoll Radiotranslucent) (Adult)	31.95 PR	Pre Con: 32.95 PR
PadPro™ Defib Pads (Physio Radiotranslucent) (Adult)	32.95 PR	Quik: 42.85 PR
PadPro™ Defib Pads (Physio Radiotranslucent) (Adult) (Pre-Con)	33.95 PR	
PadPro™ Defib Pads (Zoll Radiotransparent) (Adult)	30.50 PR	Pre Con: 31.50 PR
PadPro™ Defib Pads (Physio Radiotransparent) (Adult)	34.95 PR	Pre Con: 35.95 PR
Kendall Defib Electrodes (Physio, MediTrace Fast Patch Style)	24.95 PR	241.95 (10PR/BX)
Kendall Defib Electrodes (ZOLL MediTrace Stat Padz Style)	46.05 PR	446.65 (10PR/BX)
Kendall Defib Pads (Physio, Quick Combo w/ 40" Cable)	38.79 PR	376.25 (10PR/BX)
Kendall Defib Pads (Philips/Laerdal FR2)	36.25 PR	351.60 (10PR/BX)
Kendall Defib Pads (Hewlett Packard Code Master)	22.75 PR	220.65 (10PR/BX)
Curaplex Defib Pads (Leads Out Direct)(Zoll, Physio, Philips)	25.55 PR	
Curaplex Defib Pads (Multifunction) (Physio, Zoll, Philips)	Adult: 24.55 PR	Ped.: 24.55 PR
Curaplex Defib Pads (Multifunction) (Cardiac Sci., Powerheart G3)	29.45 EA	
Cardiac Science Adult Defib Pads (For Powerheart AEDS)	41.95 PR	Polar: 49.95 PR
Cardiac Science Pediatric Defib Pads (For Powerheart AEDS)	99.95 PR	
Conmed Defib Pads (Zoll Direct Connect R2 Multifunction) (Adult & Ped.)	32.50 PR	315.25 (10 PR/BX)
Conmed Defib Pads (Physio Direct Connect R2 Multifunction) (Adult & Ped.)	34.49 PR	334.50 (10PR/BX)
Conmed Defib Pads (Philips Direct Connect R2 Multifunction) (Adult & Ped.)	Barrel: 27.89 PR	270.50 (10PR/BX)
Conmed Defib Pads (Philips Direct Connect R2 Multifunction) (Adult & Ped.)	Plug: 31.95 PR	309.95 (10PR/BX)
Conmed Defib Pads (Welch Direct Connect R2 Multifunction) (Adult & Ped.)	26.25 PR	254.60 (10PR/BX)
Philips Heartstart FR2 Defib Pads (Adult)	42.89 PR	172.95 (5PR/BX)
Philips Heartstart FR2 Defib Pads (Pediatric)	108.75 PR	
Philips Heartstart SMART Pads Cartridge for Onsite	Adult: 62.75 PR	Ped.: 96.45 PR
Philips Heartstart MRx Multifunction Pads (Barrel Style)	Adult: 26.65 PR	Ped.: 27.95
Philips Heartstart Barrel Style Cable (Use w/ Barrel Style Pads)	138.95 EA	
Philips Heartstart MRx Multifunction Pads (Plug Style)	Adult: 31.95 PR	Ped.: 32.95 PR
Philips Heartstart FRx Defib Pads (Adult)	54.25 PR	
Philips Heartstart FR3 Smart Pads III	45.95 PR	