

PATIENT TRANSPORTATION EQUIPMENT

(INCLUDING, BUT NOT LIMITED TO):

STAIR CHAIRS & WHEEL CHAIRS

STRETCHERS

SOFT STRETCHERS PATIENT

SHIFTERS

COTS & ACCESSORIES

PEDIATRIC TRANSPORT

PATIENT RESTRAINTS

HANDLING & COMFORT

BLANKETS, SHEETS, & PILLOWS

TRIAGE SUPPLIES

STRATODYNEINC.COM

STAIR CHAIRS & WHEELCHAIRS

EZ-Glide® Stair Chairs



Our current offering of Ferno's EZ-Glide® Stair Chairs includes the Powertraxx™, Tracked, & Evacuation models. Each model features superior maneuverability with a 500 lb. maximum load. These products are equipped with easy to use operations to assure complete user control.

Ferno® Model 40 & Model 42 Stair Chairs



Both Stair Chair models include easy maneuverability in confined areas. Model 40 includes telescoping handles for easy EMT use without interfering with patient. Model 42 features glide wheels, non-slip grip handles, soft vinyl coated cover, locking safety latch, and quick release patient restraints.

MedSource Stair Chair



The MedSource Stair Chair aides in safely removing a patient from cramped quarters. Easy to maneuver by two EMT's, the unit comes with 5" rear wheels, telescoping handles, and quick-release straps. The vinyl-coated nylon cover is easy to clean and resists contamination and stains. (350 lb. Load Capacity)

PLEASE CONTACT FOR OPTIONS & PRICING

Model 40 Stair Chair (Burgundy or Orange)	876.95 EA
Model 42 Stair Chair (Burgundy or Orange)	876.95 EA

MedSource Stair Chair (36"H x 20"W x 27"D)	541.95 EA
--	-----------

Junkin Stair & Evacuation Chairs



Junkin's Stair & Evacuation Chairs were designed for confined areas. Each model holds up to 500 lbs., features a resistant and durable vinyl cover, and secure patient straps. Options include: extended handles, 4 wheels, or both.

Junkin Stair Chair	698.00 EA
Junkin Stair Chair w/ Extended Handles	808.50 EA
Junkin Stair Chair w/ 4 Wheels	768.00 EA
Junkin Stair Chair w/ Wheels & Ext. Handles	878.00 EA

Junkin Confined Space Stair Chair



Junkin's Confined Space Stair Chair features 2" standard automotive-style patient restraint straps included with each chair. Green vinyl covers are made of 18 oz. material and are fungus and rot resistant for long-lasting durability and easy cleanup. A storage/carrying case is also available. 800 lb. Load Capacity.

Junkin Confined Space Stair Chair	725.00 EA
-----------------------------------	-----------

Junkin CON Onboard Stair Chair



Junkin's CON Onboard Stair Chairs are designed for the airline industry to transport in narrow aisles. Currently used by top Airlines in USA for non-ambulatory passengers. 800lb. Load Capacity. (34.5"H x 14.75"W x 29.5" x 9")

Junkin CON Onboard Stair Chair	857.00 EA
--------------------------------	-----------

Briggs 18" Wheelchairs



Briggs 18" Wheelchairs feature chip resistant, easy-to-clean, durable, chrome-plated, and carbon-steel frame construction. Include single axle, 24" rear mag-style composite wheels with an easy to reach "push to lock" brake. They also feature a padded armrests and embossed vinyl upholstery, swing-away footrests, and lightweight foot plates. 250 lb. capacity

Wheelchair w/ Removable Arms & Swing Footrests	323.19 EA
--	-----------

Graham-Field Aluminum Transport Chair



Designed as an affordable patient transfer device. Compact size and light weight allow easy chair storage. Features flame retardant upholstery, carry straps, adjustable footrests, full length padded arms. 17" & 19" seat widths available. Meets the requirements of Medicare code E1038. 250 lb. Load Capacity.

Aluminum Transport Chair (Red, Black, Blue)	229.95 EA
---	-----------

Traveler® SE Wheelchairs



GF® Health Products has designed the Traveler® SE Wheelchair with dual rear axles and adjustable coasters. This wheelchair can convert from a K0001 Medicare code. It features industry compatible front rigging and has vinyl upholstery. It also comes with a detachable arm and swing-away footrest.

Traveler SE® (18" x 16")	229.95 EA
Traveler SE® (20" x 16")	249.95 EA



JUNKIN EASY-FOLD WHEELED STRETCHER

Made of 1" sturdy anodized aluminum tubing frame with flush folding hinges. Four inch wheels and posts fold easily out of the way for compact storage. Cover is fabricated from tough 18 oz. green vinyl coated nylon. Cover is fungus and rot resistant and easily wiped clean. 2" wide patient restraint straps with automotive seatbelt type buckles are standard equipment.

Dimensions Open: 74"L x 20.5"W x 9"H
Dimensions Folded: 36.5"L x 20.5"W x 3"H

Junkin Easy-Fold Wheeled Stretcher

326.25 EA



JUNKIN EASY-FOLD PLAIN STRETCHER

Folds for compact storage. 2" wide patient restraint straps with standard automotive seatbelt type buckles. Easy-Fold Plain Stretcher has the same frame and cover as Wheeled Model but without wheels and posts.

Dimensions Open: 74"L x 1.25"H
Dimensions Folded: 36.5"L x 2.75"H

Junkin Easy-Fold Plain Stretcher

219.39 EA



JUNKIN EASY-FOLD WHEELED STRETCHER WITH ADJUSTABLE BACK REST

Same as the Wheeled Model with the addition of an adjustable back rest with five positions and four-wheel swivel. Also available with two posts and two wheels. Folds for compact storage. 2" wide patient restraint straps with standard automotive type buckles.

Dimensions Open: 74"L x 20.5"W x 9"H
Dimensions Folded: 36.5"L x 20.5"W x 3"H

Easy-Fold 2 Wheeled Stretcher with Adjustable Back Rest
 Easy-Fold 4 Wheeled Stretcher with Adjustable Back Rest

533.35 EA
 697.45 EA



FERNO 107 COMBINATION STRETCHER CHAIR

The Ferno Model 107 is designed to serve as a wheeled chair, a stair chair, and a flat stretcher. Sturdy, lightweight aluminum frame provides long-term durability, and the vinyl-coated nylon cover is strong and easy to clean. The Model 107 features folding head-end posts and folding foot-end wheels. Includes patient restraints.

Ferno 107 Combination Stretcher Chair (Burgundy)
 Ferno 107 Combination Stretcher Chair (Orange)

1,389.95 EA
 1,389.95 EA



FERNO 107-B4 COMBINATION STRETCHER CHAIR

The Ferno Model 107-B4 is designed to serve as a wheeled chair, a stair chair, and a flat stretcher. Sturdy, lightweight aluminum frame provides long-term durability, and the vinyl-coated nylon cover is strong and easy to clean. The Model 107-B4 features fixed head-end posts, fixed foot-end wheels, and folding foot-end handles. Includes patient restraints.

Ferno 107-B4 Combination Stretcher Chair (Orange)

1,299.95 EA



FARATEC WHEELED LITTER CARRIER

The Litter Carrier is a light, stable device, equipped with two large non-pneumatic wheels designed for transporting injured persons on rugged terrain. These flat-free tires keep the unit mobile at all times, without the need for a pneumatic pump. When the supporting legs are in the transport position the patient is easily moved by one or two persons. When the supporting legs are in the "static" position the litter carrier is a steady and stable support for the patient. When not in use, the Litter Carrier is easily folded to a compact size for storage. The EL3000 Military Litter Carrier has proved very useful for transporting and moving injured patients (or materials) inside field hospitals and field hospital areas, AMPs, or to provide access to helicopter pads or airlift areas in typical large-scale emergency or natural disaster scenarios.

Wheeled Litter Carrier (Orange)
 Wheeled Litter Carrier (Black)
 NATO Wheeled Litter Carrier (Military Olive Drab)

1,995.95 EA
 1,995.95 EA
 1,995.95 EA

SOFT STRETCHERS

Stratodyne supplies a vast assortment of Soft Stretchers. From reusable to disposable options, all of our Soft Stretchers are durable, compact, and cost effective. Some features of our different models include: multi-positioned handles for added safety and comfort, pull straps for situations where pulling the patient is safer than lifting, and compact sized models that are ideal for confined situations. No matter your needs we are certain to have a product that can meet your operation's needs.



The Perry Pouch™



Simulaids Soft Stretcher



Sit™ Stretcher



Morrison Jumbo Soft Stretcher

SOFT STRETCHERS	DIMENSIONS (LxW)	CAPACITY	PRICE
MegaMover® Portable Transport Unit	40" x 80"	1,000lbs.	25.35 EA
MegaMover® Plus Portable Transport Unit	40" x 80"	1,000lbs	33.75 EA
MegaMover® Portable Transport Unit (w/ Powergrips)	40" x 80"	1,000lbs	32.25 EA
MegaMover® Select	40" x 80"	1,000lbs.	37.95 EA
MegaMover® Transport Chair	35" x 35"	350lbs.	18.75 EA
Morrison Jumbo Soft Stretcher	60" x 80"	2,000lbs.	269.95 EA
Morrison Medium Soft Stretcher	40" x 71"	2,000lbs.	151.95 EA
Morrison Small Soft Stretcher	30" x 30"	2,000lbs.	103.50 EA
Simulaids Soft Stretcher	21" x 35"	200lbs.	80.00 EA
The Perry Pouch™	28" x 31"	500lbs.	56.95 EA
Sit™ Stretcher	20" x 34"	500lbs.	49.25 EA

PLEASE CONTACT A STRATODYNE REPRESENTATIVE TODAY FOR POTENTIAL QUANTITY DISCOUNTS!!




MegaMover® Select



MegaMover®



MegaMover® Transport Chair



MegaMover® Plus

PATIENT SHIFTERS

MedSource Total Lifter™



Total Lifter™ patient transport pad is a breakthrough device that provides an easier method for lifting, moving, and extracting victims. Ideal for rescue, medics and hospitals. The Total Lifter™ holds up to 500lbs., lifts easily with 4-8 people, features anti-tear handles, provides warmth while reducing trauma impact, minimizes infections due to anti-microbial polyester, and weights only 2lbs.

Total Lifter™ 50.25 EA

Schure Slide XLA



The Schure Slide XLA is an ultra-smooth flexible sheet used to transfer patients safely and easily from bed to stretcher and/or from stretcher to bed. The sliding action eliminates the need to lift the patient, removing the risk of possible back injury to emergency personnel. Its unique patented design allows patient transfer in a sitting or supine position.

Schure Slide XLA, 33" x 18" x 1/16"	50.25 EA
Schure Slide XLA 33" x 18" x 3/16"	84.95 EA
Schure Slide XLA 72" x 18" x 3/16"	130.79 EA

Slide-Mat Patient Transfer Device



Designed and patented patient transfer device for hospital departments, EMS, and home care. Slide-Mat is an open ended design, utilizing an endless conveyor principle, eliminating the need for lifting or log-rolling to effectively transfer or reposition any patient. Slide-Mat's slick inner surface forms a unique roller. Slide-Mat's outer surface resists staining and may be washed normally.

Slide-Mat 72"x 24"	149.95 EA
Slide-Mat 18"x 24" (Wheel Chair Model)	49.95 EA

35X PROFLEXX® COT

The Ferno 35X PROFlexx® X-Frame Ambulance Cot provides industry leading stability and safety and reduces the chance of injury to you or your patients. A rigid C-channel and carefully designed frame geometry reduce your chance of tipping or flipping the cot. And the exclusive EZ-Pull™ release handle reduces wrist strain and helps you position yourself before you assume the weight of the cot.



35 PROFlexx®	5,824.75 EA
35 PROFlexx® Cot w/ Side	6,759.95 EA
35 PROFlexx® Cot w/ Stat Trac	5,788.75 EA

35-A MOBILE TRANSPORTER™

The Ferno 35-A Mobile Transporter™ X-Frame Ambulance Cot provides exceptional value for X-Frame users. The 35-A weighs only 74lbs. and has a record of dependability with over 20 years of use by EMS services all over the world. A full range of options and accessories enhance patient care and help reduce situation that may lead to patient or operator injury. The 35-A Mobile Transporter™ includes a burgundy mattress and restraints.



Ferno 35-A Mobile Transporter™	4,137.95 EA
--------------------------------	-------------

POWERFLEXX+® POWERED COT

The POWERFlexx+® Powered Ambulance Cot delivers 700 lbs. unassisted lift capacity. The cot lifts the weight, not the operator, reducing situations where operators may injure their back. The POWERFlexx+® Powered Ambulance Cot also features advanced power management features and in-the-fastener charging to keep you powered on every call, without having to carry around extra batteries or take them out to charge.



POWERFlexx+® Cot	12,130.00 EA
POWERFlexx+® Cot w/ Rotate Down Arms	12,934.95 EA

COT SHOULDER HARNESES & STORAGE NETS

EMP COT SHOULDER HARNESS


The nylon shoulder harness restraint system incorporates two shoulder straps and a two-piece torso strap that attaches to the cot to minimize longitudinal and transverse movement. Quick press release buckles allow for trouble free application. All four straps of the harness are individually adjustable to provide proper positioning.



Economy Shoulder Strap (Burgundy or Black)	65.39 EA
Impervious Shoulder Strap (Black)	45.95 EA

HARTWELL MEDICAL CATCH-ALL™

The CATCH-ALL fits most popular cots and provides a place to store commonly used equipment in a convenient location. It does not interfere with the cot's function, but it does allow the attendant easy access to the equipment at the head of the cot. Includes a zippered pouch for secure storage.



CATCH-ALL For Model 35-A	87.85 EA
CATCH-ALL For Model 93ES Cot (19" x 18.5")	87.85 EA

EMP DELUXE STORAGE NET

The deluxe storage net provides individual zippered compartments for more detailed storage of smaller supplies and equipment. A unique dedicated divided pouch holds and dispenses your gloves by size. Color coded pouches and clear windows make finding your supplies easy. Fits Ferno model 35-A and 35-A+ cots only.



Deluxe Storage Net	49.95 EA
--------------------	----------

STORAGE NET & BASE NETS FOR STRYKER

Stryker Rugged MX-Pro StoNet	65.95 EA
Stryker Rugged MX-Pro BaseNet	58.95 EA
Stryker Rugged EZ-Pro2 StoNet	65.95 EA
Stryker Rugged EZ-Pro2 BaseNet	58.95 EA





FERNO PEDI-MATE

The PEDI-MATE® is an innovative restraint system for pediatric transportation.

Features include: Use of three restraint straps to attach quickly, easily and securely to most Ferno cots; Fully-adjusting five-point harness system fits each child firmly and provides safe control during transport; Securely holds children from 10-40 lb. (4.5-18 kg); Rolls compactly allowing for efficient storage and retrieval; Made from nontoxic, easy-to-clean vinyl.

Ferno PEDI-MATE

306.95 EA



FERNO PEDI-MATE PLUS

The Ferno PEDI-MATE Plus includes the same features as the PEDI-MATE® Restraint System and accommodates larger patients ranging in size from 10-80 lbs. uses three restraint straps to attach quickly, easily and securely to most Ferno cots.

Features Include: A Fully adjustable, five-point harness system that securely holds patients, providing safe restraint for transport; Rolls compactly for convenient storage and easy deployment; Vinyl construction is nontoxic and easy-to-clean; Three restraint straps that easily attach to any cot.

Ferno PEDI-MATE Plus

313.00 EA



FERNO NEO-MATE

The Neo Mate Pediatric Restraint System adapts any ambulance cot to safely transport children ranging in size from 7-14 lb (3-6 kg).

Features include: "Halo" Pad that keeps infant's head stable; Fully-adjustable, five-point harness system fits children firmly, and provides safe restraint during transport; Three restraint straps that easily attach to any cot; Rolls up compactly for efficient storage and deployment; Designed for Infants 3 to 4.5 kg

Ferno Neo-Mate Child Transport

313.00 EA



FERNO 676 MEDKIDS BABY BOARD

The Medkids Baby Board helps maintain cervical neutral alignment in infants with the self-contained Pneumatic Positioning Device (PPD). Accommodates the wide range of thoracic elevation requirements without extra padding and does not interfere with invasive and noninvasive equipment. Made of durable DuPont Cordura and is X-ray translucent.

Medkids Baby Board

615.95 EA



SAFEGUARD PEDIATRIC TRANSPORT

One operation tightens harness and adjusts height. Integrates easily into ambulance operations, stowing efficiently in space-limited storage areas. Color coded webbing with snap hooks make attachment to cots quick and easy. Moves with the cot backrest to accommodate multiple positions from horizontal to 70 degrees, allowing easy access to the patient during transport. Extensive, rigorous crash testing assures the product is dependable and effective for transporting children safely in ambulances.

SafeGuard Pediatric Transport

506.15 EA



EP&R AIR FILLED CHILD TRANSPORT SEAT

The revolutionary Rescu-Air™ Air Filled Child Transport Seat incorporates several key improvements over its predecessor. Such as the new Cot Harness System, enabling safe and secure anchoring of the seat front and back. As well as the new release valve, ensuring proper inflation and PSI regardless of atmospheric conditions. The Rescu-Air™ was designed for ambulances and has been used by emergency medical technicians for over 12 years. Our air filled child restraint seats are unsurpassed for safety, rapid deployment, quick installation, and convenient storage.

Rescu-Air™ Air Filled Child Transport Seat, Complete

665.95 EA

Air Pump

64.95 EA

Inflation Hose Replacement


11.25 EA

Cot Harness System

79.95 EA

PATIENT TRANSPORT EQUIPMENT

PATIENT RESTRAINTS



Dick Medical Restraints



Posey Restraints



Disposable Posey Restraints



Humane Restraints

Our supply of universal restraints secure to the side rail or mainframe of the ambulance cot to hold the patients ankles and wrists securely. Our products are made to secure the patient firmly while still providing added comfort. Various construction material options are available including leather or nylon. Also, different locking mechanisms are available depending on your agency's needs.

Dick Medical Wrist & Ankle Restraints, Pair	4.55 PR
Dick Medical Wrist & Ankle Restraints with D Ring, Pair	4.85 PR
Dick Medical Nylon Wrist Restraints, Pair	14.15 PR
Dick Medical Nylon Ankle Restraints, Pair	22.99 PR
Posey Wrist Restraint, Set of 2	30.45 EA
Posey Ankle Restraint, Set of 2	30.19 EA
Posey Wrist & Ankle, Set of 4	66.95 EA
Disposable Posey Restraints Foam Limb Holder, Pair	7.09 PR
Humane Leather Wristlet, 3"W x 12"L	84.35 PR
Humane Leather Wristlet with Foam Liner	87.29 PR
Humane Leather Ankle, 3"W x 14"L	85.29 PR
Humane Leather Ankle with Foam Liner	88.19 PR
Humane Leather Belt, 1 1/4"W x 18"L	30.00 EA
Humane Leather Belt, 1 1/4"W x 24"L	32.85 EA
Humane Leather Belt with Key Lock, 1 1/4"W x 18"L	36.79 EA

PATIENT COMFORT PRODUCTS



QT Disposable Patient Transport



Ready-Heat™ II Blankets




Cot Quilt

Stratodyne offers a variety of products for Patient Comfort. We carry many diverse lines of specialty blankets and quilts to keep patients warm and protected from the environment. In addition to blankets, we also offer disposable & reusable patient transportation methods for all patient sizes. These include transfer sheets as well as Stair Chairs. Our patient comfort line also includes a Abdominal Restraint to allow for patient stability, comfort, and safety, while at the same time, making transport easier for paramedics.

QT Disposable Patient Transport (10 Per Case)	27.95 EA	271.10 CS
Ready-Heat II, 12 Panel Blanket 34" x 48"	23.99 EA	186.15 (8/CS)
Ready-Heat II, Panel Blanket 16" x 18"	10.69 EA	127.95 (12/CS)
Patient Cover Quilt, 72"W x 96"L	99.15 EA	
Patient Cover Quilt /w Window Cutouts for Stryker Cot	112.59 EA	
Cot Quilt, Black/Blue	124.95 EA	
Cot Quilt, Burgundy/Blue	124.95 EA	
B.E.A.R. - Bariatric Equalizing Abdominal Restraint	389.95 EA	
B.E.A.R. Stair Chair	545.00 EA	

CONTACT STRATODYNE TODAY FOR POTENTIAL VOLUME DISCOUNTS



B.E.A.R. Stair Chair



Patient Cover Quilt



B.E.A.R. - Bariatric Equalizing Abdominal Restraint



Flat Sheet



Wool Non-Disposable



Polyester



Water Repellent Disposable

PATIENT BLANKETS


Stratodyne carries a huge inventory of Patient Blankets. From washable to disposable, we offer blankets of all varying dimensions depending on your needs. Blanket material options include: Polyester, Wool, Radio Translucent Composite, Mylar, Poly-Tissue, Foam Insulated, Fleece, and combination products. We also carry weather resistant and waterproof options.

VISIT STRATODYNEINC.COM FOR A COMPLETE LISTING OF PRODUCTS & SPECIFICATIONS

Curaplex Flat Sheet	50.95 EA	
Curaplex Non-Disposable Blanket	16.50 EA	399.95 (25/CS)
Curaplex Polyester Blanket (72"x 90") (Blue)	8.50 EA	102.00 (12/CS)
Curaplex Polyester Blanket (60"x 90") (Gray)	7.19 EA	57.50 (12/CS)
Graham Medical Comfort1® Blanket, Polyester (50"x 84") (Gray)	8.75 EA	84.85 (10/CS)
Graham Medical Comfort1® Blanket, Polyester (50"x 84") (White)	9.75 EA	94.55 (10/CS)
Thermo-Lite® Infant Papoose	4.79 EA	
Thermo-Lite® Transport Cocoon™, Pediatric	12.60 EA	
Thermo-Lite® Transport Cocoon™, Adult	16.95 EA	
EMP Emergency Blanket (84"x 52")	2.15 EA	
Curaplex Mylar Blanket (84"x 52")	1.10 EA	
Curaplex Yellow Blanket 90"x 56"	3.95 EA	
EMP Disposable Blanket (90"x 58")	5.99 EA	103.15 (18/CS)
Curaplex Reusable Fleece Blanket 90"x60"	13.79 EA	137.90 (10/CS)
EMP Heavy Weight Wool Blanket (80"x 60")	18.89 EA	
EMP Disposable Polyester Blanket (80"x 60")	7.39 EA	179.25 (25/CS)



Curaplex Cot Sheet



Value Bedding Pack Plus



Snugfit Expandable Mattress Cover



FlexDrape™ Disposable

COT SHEETS


We currently offer a vast assortment of styles, fabrics, and sizes of Cot Sheets and Cot Sheet Kits. Our sizing options (including Snug fit expandable sheets) are certain to offer an option for any size or brand of cot your agency currently utilizes. Our inventory also includes sheeting for stretchers, patient privacy sheets, and utility sheets. Disposable and reusable options are available. Our bedding packs include (1) fitted sheet, (1) flat sheet, (1) protective underpad, and (1) pillow case. Our kits are also available in disposable and reusable options.

CONTACT A STRATODYNE REPRESENTATIVE TO DISCUSS POTENTIAL VOLUME DISCOUNT PRICING FOR THESE PRODUCTS.

MedSource Fitted Cot Sheet, Elastic Ends		64.95 (50/CS)
MedSource Fitted Cot Sheet for Stryker XPS		86.65 (30/CS)
Curaplex Cot Sheet 73"x 30" (Elastic Ends)		49.95 (50/CS)
Curaplex Cot Sheet 73"x 30" (Fitted Square Ends)		45.25 (50/CS)
Curaplex Cot Sheet 73"x 30" (Elastic Ends) (Ultra Fluid Resistant)		126.50 (50/CS)
Curaplex Cot Sheet 73"x 30" (Elastic Ends) (Heavy Duty)		79.50 (50/CS)
EMP Disposable Cot Sheet 90"x 40"		31.95 (50/CS)
FlexDrape™, 50" x 84"		74.59 (50/CS)
FlexDrape™, 35" x 50"		35.95 (50/CS)
Graham MegaMover 350i Barrier Sheet 84"x 40"		62.95 (50/CS)
Patient Privacy Sheets 84"x 40"		62.95 (50/CS)
EMP Value Bedding Pack	3.90 PK	94.55 (25/CS)
Graham Standard InstaKit	4.05 PK	99.95 (24/CS)
EMP Disposable Bedding Pack	3.49 EA	84.60 (25/CS)
EMP Disposable Bedding Pack for Power Cots	4.49 EA	108.88 (25/CS)
MedSource Linen Set	3.80 EA	92.19 (25/CS)
Medsource Linen Set (Heavy Duty)	4.25 EA	102.95 (25/CS)
SungFit® EMS Barrier Fitted Sheet* (All Sizes)		Pricing Between 64.06 - 91.89 (50/CS)
PowerFit® Sheet* Nonwoven (Blue) 89"x 33"		91.50 (50/CS)
PowerFit® Barrier Sheet* Nonwoven/Poly (Blue) 89"x 33"		107.95 (50/CS)
Snugfit Expandable Mattress Cover, Blue		73.50 (30/CS)
EMP Fitted Cot Sheets 72"x 30"		47.69 (50/CS)
Disposable Fluid Resistant Fitted Cot Sheets 72"x 30"		66.45 (50/CS)
Disposable Extra Heavy Fitted Stretcher Sheets		96.95 (50/CS)




FlexAir™ Disposable Pillows



Disposable Pillows



InhibX® Antimicrobial Pillowcase



EMP Disposable Pillow Cases

COT PILLOW & PILLOW CASES

Currently stocking a variety of reusable & disposable Cot Pillows and Pillow Cases. Versatile, comfortable, non-allergenic, tough, and cost effective options are available to meet your specific needs.

FlexAir Disposable Pillow , 14.5" x 10.5"		34.99 (50/CS)
FlexAir Disposable Pillow, 19" x 12.5"		41.60 (50/CS)
E-Space Disposable Pillows 17"x 23"	5.40 EA	64.80 (12/CS)
Morrison Staph-Chek Reusable Pillows 16"x 21"	11.09 EA	64.59 (6/CS)
InhibX® Pillowcases 21"x 30"		43.25 (80/CS)
EMP Disposable Pillow Cases 21"x 30"		31.95 (100/CS)
EMP Disposable Pillow Cases 21"x 30" (Fluid Resistant)		41.95 (100/CS)

TRIAGE SUPPLIES



PATIENT MEDICATION BAG

This handy bag solves the problem of where to put your patients medicines on the way to the ER. Constructed of high grade vinyl with a seal able flap. The bag front is white and the back is clear so providers are able to see the names of medications that are sealed in pouch.

Patient Medical Bag (8.5"x 9")	1.39 EA
--------------------------------	---------



TRIAGE KIT COMPLETE (SMALL)

(1) Commander Vest, (1) Triage Vest, (1) Treatment Vest, (1) Transport Vest, (1) Staging Vest. (1 Posse BX Plain, (25) Triage Bags, (1) Triage Tape of each: Red, Yellow, Green, Black/White, (1) San Marcos Cert Red Bag.

Triage Kit Complete (Small)	263.95 EA
-----------------------------	-----------



THOMAS TRIAGE PACK

Trauma shears, N95 Mask (1) Large, (1) Medium, Antiseptic Cleaning Towelettes, 2" Cloth Tape, Penlight, Sharpie Pen. Also includes: (6) Rolls of Flagging (Red, Yellow, Black, Blue, Green, & Orange). Additional triage and treatment items can be added to this pack per request.

Thomas Triage Pack	85.75 EA
--------------------	----------



STANDARD MEDICAL EMER. FIELD TRIAGE TAG

This tag is used extensively by U.S. military units, as well as numerous Federal, state and local response agencies, including the Dept. of Defense, Dept. of Corrections, state and local Emer. Mgmt. Agencies and Operations Centers, etc.

Standard Field Triage Tags	55.45 (50/PK)
California Fire Chiefs Triage Tags	81.95 (50/PK)



BODY BAGS

Body Bag measures 36"x 90". Pediatric measures 33"x40. Premiums measures 85"x 36". Disaster pouch measures 36"x 90". Weight limit is 300lbs.

Body Bag with Toe Tags, 8-10 mm Thick	13.49 EA
Pediatric Post-Mortem Kit	10.99 EA
Premium Body Bag	20.85 EA
Disaster Pouch, 16-20 mm Thick	41.95 EA



ACCOUNTABILITY/COMMAND BOARD

This board helps keep track of your people and organizes your incident command structure. Complete dimensions of the board are 19"x 25"x 5") The Accountability/Command board is a must for organizing large operations .

Accountability / Command Board	179.95 EA
--------------------------------	-----------



TRIAGE FLAG KITS

Our triage flags are sewn entirely out of U.S. Oxford nylon. The flag kit comes with red, yellow and green flag/poles, three steel bases and a carrying pouch. A must for every emergency vehicle. Automatically deploys to 54" High.


TRIAGE Flag Kit (Red, Yellow, Green)	117.95 EA
TRIAGE Flag Kit (Red, Yellow, Green, Black)	146.95 EA



EASY COT

Easy to set up and easy to store full size cot. 32" cot width with a set up time of under 1 minute. Welded steel legs, multiple fabric reinforcement, assures that it is incredible strong and durable no matter your needs. Comes with storage bag as show.

Easy Cot	75.75 EA
----------	----------



TRIAGE TARP

The triage tarps are made from heavy coated nylon oxford cloth for years of use. Heavily reinforced corners with 1000d Cordura nylon, & feature powder coated brass grommets to use as tie downs, or suspension points. Color coded to match S.T.A.R.T. requirements, the tarps measure a full 10' x 12'.

Triage Tarp (Red, Yellow, Green, Black)	116.95 EA
---	-----------



DELUXE TRIAGE BELT

Deluxe Triage belt organizes your gear in a totally new way. The slim profile hugs your body, and the totally enclosed tape won't tangle or flap. The belt sports a large front pocket, with two outside pockets, and a lanyard keeper. Lifetime guaranteed.

Deluxe Triage Belt	48.95 EA
--------------------	----------



TRIAGE BELT

Unique belt designed for use within the S.T.A.R.T. triage system. Provides a organized and effective method for patient taping. This product will eliminate confusion that often occurs due to double tagging. Belt tape rolls measure 300' in length.

TRIAGE Belt	22.95 EA
Belt Tape (Red, Yellow, Green, Blue, Black)	2.59 EA



MCI KIT

Includes: (1) Command Vest, (1) Safety Vest, (1) Triage, (1)Treatment Vest, (1) Transport Vest, (1) Medical Vest. The kit also includes: a start triage belt, and has several internal pockets for all necessary triage items. Dimensions are 16"x16"x8").

MCI Kit Complete	296.95 EA
------------------	-----------

POWERFlexx[®]+

ADVANCED POWER MANAGEMENT

The NEW **POWERFlexx+** System delivers advanced power management features and **in-the-fastener charging** with the ICS Integrated Fastener Charging System, to keep you powered on every call without having to carry around extra batteries or take them out to charge.



700 LB UNASSISTED LIFTING

The **POWERFlexx+** Powered Ambulance Cot also features 700 lb unassisted lift capacity. The cot lifts the 700 lb. weight, not the operators, reducing situations where operators may injure their back.



FEATURES

POWERFUL

- ▶ The ICS - Integrated Charging System is built into the cot fastener and keeps the battery powered for an unlimited number of runs, eliminating the need to carry around extra batteries or take them out to charge
- ▶ Sleep mode preserves battery life and wakes up at the touch of the activation pad, so a strong charge is always there when you need it.
- ▶ Hydraulic Actuator Unit with dual rail system provides smooth and stable activation.
- ▶ 700 lb unassisted lifting capacity, the highest in the industry.

SAFE

- ▶ 700 lb lifting capacity reduces situations where operators may injure their back.
- ▶ No pinch points anywhere along the mainframe of the cot.
- ▶ Wide Wheel Base, Placement of Motor & Battery, and proper X-frame positioning create lower center of gravity and greater stability.

MANEUVERABLE

- ▶ Five Position Drop-frame and multi-position Foot-end Lift Bar provide better leverage for lifting and loading and accommodates movement in confined spaces.
- ▶ Infinite bed height adjustment from 12" to 41" giving you the greatest range and flexibility to transfer the patient on to the cot.

COMFORTABLE

- ▶ Improved Patient Surface increases surface length by 6¼ inches to accommodate taller patients and keep the patient's feet away from operators.
- ▶ Dual Lifting Cylinders and Rails provide smooth and even lifting for greater patient sense of security.

FLEXIBLE BARIATRICS

- ▶ No need to have a dedicated, bariatric only, cot. The **POWERFlexx+** features 1000 lb capacity and full 360° maneuverability in the lowest position and quickly converts to a full bariatric cot with accessories such as the removable LBS Jr. Bariatric Board, Manta Transfer Sheet, Push Pull Handles, and Tow Package.

SPECIFICATIONS

Length:	Maximum:	83" (211 cm)
	Minimum:	61" (155 cm)
Width:		24" (61 cm)
Max Loading Height:		34" (84 cm)
Weight:		132 lb (60 kg)
Lifting Capacity:		700 lb (318 kg)
Load Limit:	Loading Position:	700 lb (318 kg)
	Lowest Position:	1000 lb (450 kg)

VISIT STRATODYNEINC.COM FOR MORE INFORMATION & SIMILIAR PRODUCTS